## FORMULATION HANDBOOK



# 4282-055 **Key Words:** Diclofenac Sodium, Good

Flowability, Direct Compression

JRS Products: PROSOLV® EASYtab SP

## **Diclofenac Sodium Direct Compression**

Page 1/2

## **Summary**

Diclofenac is a is a non-steroidal anti-inflammatory drug (NSAID). It is used as a pain-reliever and for its antipyretic (fever-reducing) effects.

Common dosages in tablets are 25 and 50 mg Diclofenac Sodium.

#### **Formulation**

	Active content [mg]	mg/tablet	Contribution [%]
Diclofenac Sodium	25.0	25.0	12.5
PROSOLV® EASYtab SP		175.0	87.5
Total		200.0	100.0

#### **Procedure**

#### **Blending:**

The API was blended with PROSOLV® EASYtab SP for 15 minutes. The powder mixture was ready for direct compression.

Equipment:	
Tablet Press:	Kilian Pressima, 8 mm punch
Turbula Mixer:	Type T2A
Hardness Tester:	Schleuniger Tablet Tester 6D
Dissolution Tester:	Pharmatest PTW 2
Spectrophotometer:	Cecil CE 1021

### **Tablet Characteristics**

Compaction force:	3.0 kN	
Tablet diameter:	8 mm	
Tablet weight:	200 mg	
Crushing strength:	70 N	
Crushing strength: Friability:	0.0 %	
Disintegration time:	20 s	



JRS PHARMA GMBH & CO. KG Excipients • CoatingsBiopharma Services • Technical Services

73494 Rosenberg (Germany) Phone: + 49 (0) 7967 / 152 312 Fax: + 49 (0) 7967 / 152 345

**WORLDWIDE HEADQUARTERS** 

ExcipientsService@jrspharma.de

# FORMULATION HANDBOOK



# 4282-055 Key Words: Diclofenac Sodium, Good

Flowability, Direct Compression

JRS Products: PROSOLV® EASYtab SP

# **Diclofenac Sodium Direct Compression**


Page 2/2

**Dissolution Test:** 

Dissolution Medium: 900 mL 0.1 N HCl, 37°C, n=6

Samples were taken after 5, 10, 15, 30, 60, 120 and 180 minutes.

The determination of the Diclofenac was done by UV spectroscopy at  $\lambda = 280$  nm.



**Diagram 1:**Typical dissolution profile of a Diclofenac Sodium tablet. Produced according to the above formulation.

**Disclaimer:** The information provided is based on thorough research and is believed to be completely reliable. Application suggestions are given to assist our customers, but are for guidance only. Circumstances in which our material is used vary and are beyond our control. Therefore, we cannot assume any responsibility for risks or liabilities, which may result from the use of this technical advice.



• Excipients • Coatings

• Biopharma Services • Technical Services

WORLDWIDE HEADQUARTERS
JRS PHARMA GMBH & CO. KG

73494 Rosenberg (Germany) Phone: + 49 (0) 7967 / 152 312 Fax: + 49 (0) 7967 / 152 345

ExcipientsService@jrspharma.de www.irs.de